

几种 AAV 在脑中的表达比较

参考文献: Analysis of Transduction Efficiency, Tropism and Axonal Transport of AAV Serotypes 1, 2, 5, 6, 8 and 9 in the Mouse Brain

AAV1, AAV2, AAV5, AAV6, AAV8, AAV9 这 6 种 AAV 是目前最常用的 AAV 类型。

Striatum (纹状体): AAV5 > AAV9 > AAV8 > AAV1 > AAV6 > AAV2;

Hippocampus (海马): AAV9 > AAV5 > AAV8 > AAV1 > AAV6 > AAV2;

Auditory cortex (听觉皮层): AAV9 > AAV8 > AAV6 > AAV5 > AAV1 > AAV2.

