

## 客户用百恩维提供的 AAV5 注射肺发的文章

文章标题: miR-29a-3p Attenuates Hypoxic Pulmonary Hypertension by Inhibiting Pulmonary Adventitial Fibroblast Activation

期刊:

Hypertension.2015;65:00-00.DOI:10.1161/HYPERTENSIONAHA.114.04600

注射方法: Rats were anesthetized with pentobarbital (30 mg / kg, i.p.), and given a double intratracheal instillation of AAV5-miR-29a or AAV5-control at day-0 and day-14 (for preventive experiment) or a single intratracheal instillation at day-28 after hypoxia (for therapeutic experiment) ( $1 \times 10^{12}$  v.g / kg / each instillation) according to the grouping.